

# CASE STUDY



Dr Annabelle Trainor BSc, DVM

CROWN EQUINE VETERINARY SERVICES – MELBOURNE

## EQUINE - Non healing leg wound 6 months



### K-LASER THERAPY

Twice a week for 8 weeks then once a week for a month until healing

Settings for chronic wounds

4 wavelengths - 660 -800 -905 -970nm

Average Power 3 watts

11 phases with

Continuous wave and

Pulsing 10 to 10.000Hz



1rst K-Laser  
Session

Healing

25th Oct 2019



15th Nov 2019



16 Dec 2019



17th Jan 2020

